

Technical Data Sheet

ACS Material Single Layer Graphene Oxide Dispersion

Table of Contents

- 1 Preparation Method
- 2 Characterizations
- 3 Application Fields

Contact Information:

Manufacturer: ACS Material, LLC. Address: 959 E Walnut St., Suite 100 Pasadena, CA 91106, USA Phone: (866)-227-0656 Fax: (781)-518-0284 E-Mail: contact@acsmaterial.com Revision: 062017

1. Preparation Method

All dispersions are made from our Single Layer Graphene Oxide (H Method).

2. Characterizations

New SKU:	GNO1W001	GNODAW01 (New Product)	GNO1W005	GNOD1E01	GNOD1W01
Concentration:	10mg/ml.	5mg/ml.	5mg/ml.	5mg/ml.	5mg/ml.
	100ml per	1L per bottle	100ml per	100ml per	100ml per
	bottle (1g)	(5g)	bottle (0.5g)	bottle (0.5g)	bottle (0.5g)
Solvent:	DI Water	DI Water	DI Water	Ethanol	DI Water
Flake size:	0.5-2.0 μm	0.5 - 2.0 μm	~ 500 nm	0.5-2.0 μm	0.5 - 2.0 μm
Thickness:	0.6-1.2 nm	0.6-1.2 nm	0.6-1.2 nm	0.6-1.2 nm	0.6-1.2 nm
Single-layer	>80%	>80%	>80%	>80%	>80%
Ratio:					
Color:	Brown/Black	Brown/Black	Brown/Black	Brown/Black	Brown/Black



Typical TEM Image of ACS Material Graphene Oxide Dispersion



AFM Analysis of ACS Material Graphene Oxide Dispersion

3. Application Fields

- 1) Catalyst
- 2) Supercapacitors
- 3) Solar energy
- 4) Graphene semiconductor chips
- 5) Conductive graphene film
- 6) Graphene computer memory
- 7) Biomaterials
- 8) Transparent conductive coatings

Disclaimer: ACS Material, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. ACS Material makes no representations or warranties either express or implied, regarding the suitability of the material for any purpose or the accuracy of the information contained within this document. Accordingly, ACS Material will not be responsible for damages resulting from use of or reliance upon this information.